

FluoroQuest™ Fluorescence Quantum Yield Determination Kit *Optimized for Bioconjugates*

Catalog number: 20

Unit size: 1 kit

Component	Storage	Amount
Component A: FRS520	Freeze (<-15 °C), Minimize light exposure	1 mL
Component B: FRS570	Freeze (<-15 °C), Minimize light exposure	1 mL
Component C: FRS610	Freeze (<-15 °C), Minimize light exposure	1 mL
Component D: FRS670	Freeze (<-15 °C), Minimize light exposure	1 mL

OVERVIEW

When a fluorophore absorbs a photon of light, an energetically excited state is formed. The fate of this species is varied, depending upon the exact nature of the fluorophore and its surroundings, but the end result is deactivation (loss of energy) and return to the ground state. The main deactivation processes which occur are fluorescence (loss of energy by emission of a photon), internal conversion and vibrational relaxation (non-radiative loss of energy as heat to the surroundings), and intersystem crossing to the triplet manifold and subsequent non-radiative deactivation. The fluorescence quantum yield is the ratio of photons absorbed to photons emitted through fluorescence. In other words the quantum yield gives the probability of the excited state being deactivated by fluorescence rather than by another, non-radiative mechanism. The kit provides all the essential components for fluorescence quantum yield determination for biological conjugates. It is optimized for determining the fluorescence quantum yields of fluorescent protein conjugates, peptides, nucleotides and nucleic acids.

AT A GLANCE

Experimental considerations

The most reliable method for recording fluorescence quantum yield is the comparative method of Williams *et al.* It involves the use of well characterized standard samples with known fluorescence quantum yield values. Essentially, solutions of the standard and test samples with identical absorbance at the same excitation wavelength can be assumed to be absorbing the same number of photons. Hence, a simple ratio of the integrated fluorescence intensities of the two solutions (recorded under identical conditions) will yield the ratio of the quantum yield values. Since the fluorescence quantum yield for the standard sample (Φ_s) is known, it is trivial to calculate the fluorescence quantum yield for the sample (Φ_x).

The standard reference should be chosen to have certain absorbance at the excitation wavelength of choice for the test sample, and, if possible, emit in a similar region to the test sample. Standard 10 mm path length fluorescence cuvettes are sufficient for running the fluorescence measurements. For small amount samples, either a microplate reader or nanodrop device can be used. The absorbance in the 10 mm cuvette should be <0.2. Above this level, non-linear effects may be observed due to inner filter effects, and the resulting quantum yield values may be inaccurate.

This kit is only designed for a one-point quick estimation of the fluorescence quantum yield of a biological conjugate. For more accurate deamination of fluorescence quantum yield there are many factors that have a significant effect on fluorescence quantum yields.

KEY PARAMETERS

Instrument: Spectrofluorometer
 Excitation: See Table 1
 Emission: See Table 1

Instrument: Fluorescence microplate reader
 Excitation: See Table 1

Emission: See Table 1
 Recommended plate: Solid black

PREPARATION OF STOCK SOLUTIONS

Unless otherwise noted, all unused stock solutions should be divided into single-use aliquots and stored at -20 °C after preparation. Avoid repeated freeze-thaw cycles.

Make 1 mM stock solutions of the selected standard and test samples using a proper solvent (such as water or DMSO).

Table 1. Kit components with respective quantum yield in water and similar wavelength dyes for reference standard selection.

Components	Fluorescence Quantum Yield in Water (Φ_s Value)	Conjugates Labeled with the Following Dyes or Dyes of Similar Wavelength (for Reference Standard Selection)
Reference A	0.98	FAM, 6-TET, 6-HEX, 6-JOE, FITC, Cy ² , Alexa Fluor [®] 488 and 514, iFluor [™] 488 and 514, DyLight [™] 488, or other dyes that have an emission of 500 ± 50 nm
Reference B	0.20	Cy ³ , Alexa Fluor [®] 514, 532, 546 and 555, iFluor [™] 514, 532 and 555, DyLight [™] 555, TRITC, or other dyes that have an emission of 550 ± 50 nm
Reference C	0.44	Texas Red [®] , Texas Red [®] -X, Alexa Fluor [®] 594, iFluor [™] 594, California Red [™] , DyLight [™] 594, or other dyes that have an emission of 600 ± 50 nm
Reference D	0.24	Cy ⁵ , Cy5.5 [®] , Cy7 [®] , Alexa Fluor [®] 633, 647, 700 and 750, iFluor [™] 633, 647, 700 and 750, DyLight [™] 650, 680 and 755, or other dyes that have an emission of 650 ± 50 nm

SAMPLE EXPERIMENTAL PROTOCOL

- Record the absorbance spectra of the selected standard and test samples (0.1 to 5 μ M) using the diluted solutions prepared from their stock solutions with an aqueous buffer (such as PBS or TRIS).
- Select the excitation wavelength based on the absorption maxima. We recommend to use this following formula to select the desired excitation wavelength:

$$\text{Desired excitation (nm)} = \text{shorter absorption maximum}^* - 50\text{nm}$$
 *Selected from either the reference standard or test sample. For example, if your test sample has absorption maximum at 500 nm, and the selected reference standard has the absorption maximum at 490 nm, you should select 440 nm (490 nm - 50 nm) as excitation wavelength.
- Adjust the concentrations of the selected reference standard and test samples to have the same absorbance (between 0.2 to 0.8) at the selected excitation wavelength.

Note It might be easier to choose the isosbestic point of the reference standard and test sample absorption spectra.

- Dilute both the selected reference standard and test sample by 10 times (must be diluted by the same proportion) using the same buffer or pure water.
- Record the fluorescence spectra of the same diluted solutions in the 10 mm fluorescence cuvette using a fluorescence spectrophotometer or a microplate reader that can generate corrected fluorescence spectrum.

Note All the fluorescence spectra must be recorded with the same instrument settings. Changing instrument settings between samples will invalidate the quantum yield measurement.

EXAMPLE DATA ANALYSIS AND FIGURES

- Calculate the integrated fluorescence intensity (the area of the fluorescence spectrum) from the fully corrected fluorescence spectrum (see your instrument's instructions).
- Calculate the fluorescence quantum yield using the following equation: $\Phi_r = \Phi_s \cdot A_r/A_s$
 Φ_r and Φ_s are the fluorescence quantum yields of the selected reference standard and test sample respectively. A_r and A_s are the integrated fluorescence intensity (the curve area) of selected reference standard and test sample respectively.




Figure 1.

Principle of fluorescence. Electrons are excited to a higher energy level by external source. Upon return to their ground state, a set quanta of photons are release proportional to the energy loss by electrons. This release of photons represents the fluorescence emission. The fluorophore's quantum yield is the ratio of its emitted photons to the photons it absorbed.

DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email info@aatbio.com if you have any questions.