

## Cell Meter™ Cell Viability Assay Kit \*Blue Fluorescence with 405 nm Excitation\*

Catalog number: 22784  
Unit size: 500 Tests

Component	Storage	Amount
Component A: CytoCalcein™ Violet 450, AM	Freeze (<-15 °C), Minimize light exposure	5 vials, lyophilized
Component B: DMSO	Freeze (<-15 °C)	1 vial (200 µL)
Component C: Assay Buffer	Freeze (<-15 °C)	1 bottle (50 mL)

### OVERVIEW

There are a variety of parameters that can be used to monitor cell viability. The proprietary violet laser -excitable fluorescent dye used in the kit is a hydrophobic compound that easily permeates intact live cells and gets enhanced fluorescence upon entering into live cells. The hydrolysis of the non-fluorescent substrate by intracellular esterases generates a strongly blue fluorescent product that is well-retained in the cell cytoplasm. The blue fluorophore generated by the esterase hydrolysis of the non-fluorescent substrate has the spectral properties of fluorescein. When excited at 405 nm, the fluorophore emits intense blue fluorescence at ~450 nm. The kit provides all the essential components with an optimized cell-labeling protocol for fluorescence microplate assays. This Cell Meter™ Cell Viability Assay Kit provides an effective tool of labeling cells for fluorescence flow cytometry, microplate and microscopic investigations of cellular functions. It is useful for a variety of studies, including cell adhesion, chemotaxis, multidrug resistance, cell viability, apoptosis and cytotoxicity. The kit is suitable for proliferating and non-proliferating cells.

### AT A GLANCE

#### Protocol summary

1. Prepare cells with test compounds
2. Add the same volume of CytoCalcein™ Violet 450, AM working solution (100 µL/well/96-well plate or 25 µL/well/384-well plate)
3. Incubate at room temperature or 37°C for 1 hour
4. Monitor fluorescence intensity (bottom read mode) at Ex/Em = 405/460 nm (Cutoff = 435 nm)

**Important** Thaw one of each kit component at room temperature before starting the experiment.

### KEY PARAMETERS

Instrument:	Fluorescence microplate reader
Excitation:	405 nm
Emission:	460 nm
Cutoff:	435 nm
Recommended plate:	Black wall/clear bottom
Instrument specification(s):	Bottom read mode

### PREPARATION OF STOCK SOLUTIONS

Unless otherwise noted, all unused stock solutions should be divided into single-use aliquots and stored at -20 °C after preparation. Avoid repeated freeze-thaw cycles.

#### 1. CytoCalcein™ Violet 450, AM stock solution:

Add 20 µL of DMSO (Component B) into the vial of CytoCalcein™ Violet 450, AM (Component A), and mix well to make CytoCalcein™ Violet 450, AM stock solution.

**Note** 20 µL of CytoCalcein™ Violet 450, AM stock solution is enough for one plate. Protect from light. For storage, seal tubes tightly.

### PREPARATION OF WORKING SOLUTION

Add the whole content (20 µL) of CytoCalcein™ Violet 450, AM stock solution into 10 mL of Assay Buffer (Component C), and mix well to make CytoCalcein™ Violet 450, AM working solution. This CytoCalcein™ Violet 450, AM working solution is stable for at least 2 hours at room temperature.

**Note** If the cells such as CHO cells contain organic-anion transporters which promote the leakage of the fluorescent dye over time, a probenecid stock solution should be prepared and added to the loading buffer at a final in-well working concentration ranging from 1 to 2.5 mM. Aliquot and store the unused probenecid stock solution at ≤ -20 °C.

### PREPARATION OF CELL SAMPLES

For guidelines on cell sample preparation, please visit <https://www.aatbio.com/resources/guides/cell-sample-preparation.html>

### SAMPLE EXPERIMENTAL PROTOCOL

1. Treat cells with test compounds as desired.

**Note** It is not necessary to wash cells before adding compounds. However, if tested compounds are serum sensitive, growth medium and serum factors can be aspirated away before adding compounds. Add 100 µL/well (96-well plate) and 25 µL/well (384-well plate) of 1X Hank's salt solution and 20 mM Hepes buffer (HHBS) or the buffer of your choice after aspiration. Alternatively, cells can be grown in a serum-free media.

2. Add 100 µL/well (96-well plate) or 25 µL/well (384-well plate) of CytoCalcein™ Violet 450, AM working solution.
3. Incubate the plate at room temperature or 37°C for 1 hour, protected from light. (The incubation time could be from 15 minutes to overnight. We got the optimal results with the incubation time less than 4 hours).

**Note** The appropriate incubation time depends on the individual cell type and cell concentration used. Optimize the incubation time for each experiment.


**Note** DO NOT wash the cells after loading.

**Note** For non-adherent cells, it is recommended to centrifuge cell plates at 800 rpm for 2 minutes with brake off after incubation.

4. Monitor the fluorescence intensity with a fluorescence plate reader (bottom read mode) at Ex/Em = 405/460 nm (Cutoff = 435 nm).

### EXAMPLE DATA ANALYSIS AND FIGURES

The reading (RFU) obtained from the blank standard well is used as a negative control. Subtract this value from the other standards' readings to obtain the baseline corrected values. Then, plot the standards' readings to obtain a standard curve and equation. This equation can be used to calculate CHO-K1 Cells samples. We recommend using the Online Linear Regression Calculator which can be found at:



**Figure 1.** CHO-K1 cell number response was measured with Cell Meter™ Cell Viability Assay Kit. CHO-K1 cells at 0 to 5,000 cells/well/100  $\mu$ L were seeded overnight in a Costar black wall/clear bottom 96-well plate. The cells were incubated with 100  $\mu$ L/well of CytoCalcein™ Violet 450, AM dye-working solution for 1 hour at room temperature. The fluorescence intensity was measured at Ex/Em = 405/460 nm with NOVOstar instrument (from BMG Labtech). The fluorescence intensity was linear ( $R^2 = 1$ ) to the cell number as indicated.

#### DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.