

Cell Meter™ Phosphatidylserine Apoptosis Assay Kit *Blue Fluorescence Excited at 405 nm*

Catalog number: 22835

Unit size: 100 Tests

Component	Storage	Amount
Component A: Apopxin™ Violet 450 (100X stock solution)	Refrigerate (2-8 °C), Minimize light exposure	1 vial (200 µL)
Component B: Assay Buffer (4 °C)	Refrigerate (2-8 °C)	50 mL
Component C: 100X Propidium Iodide	Freeze (<-15 °C), Minimize light exposure	1 vial (100 µL)

OVERVIEW

Our Cell Meter™ assay kits are a set of tools for monitoring cell viability. There are a variety of parameters that can be used for monitoring cell viability. This particular kit is designed to monitor cell apoptosis through measuring the translocation of phosphatidylserine (PS). In apoptosis, PS is transferred to the outer leaflet of the plasma membrane. The appearance of phosphatidylserine on the cell surface is a universal indicator of the initial/intermediate stages of cell apoptosis and can be detected before morphological changes can be observed. Our proprietary Apopxin™ PS sensor used in this kit is small molecule-based PS sensor. The PS sensor has blue fluorescence upon binding to membrane PS. The PS stain used in the kit has the spectral properties similar to those of Pacific Blue® (Pacific Blue® is the trademark of Invitrogen). The blue fluorescent stain is well excited with the Violet Laser at 405 nm, and emits intense blue fluorescence at ~450 nm. The kit is optimized to be used with a flow cytometer equipped with a Violet Laser. It is particularly suitable for multicolor flow cytometric analysis of cells. This kit has been used for flow cytometric analysis of cells in multicolor applications in combination with fluorescent antibodies.

AT A GLANCE

Protocol summary

1. Prepare cells with test compounds (200 µL/sample)
2. Add Apopxin™ Violet 450 assay solution
3. Incubate at room temperature for 30 - 60 minutes
4. Analyze cells using flow cytometer with 450/40 nm filter (Pacific Blue channel) or fluorescence microscope with Violet filter

Important Thaw 100X Propidium Iodide (Component C) at room temperature before starting the experiment.

KEY PARAMETERS

Instrument: Flow cytometer
 Excitation: 405 nm laser
 Emission: 450/40 nm filter
 Instrument specification(s): Pacific Blue channel

Instrument: Fluorescence microscope
 Excitation: Violet filter
 Emission: Violet filter
 Recommended plate: Black wall/clear bottom

PREPARATION OF CELL SAMPLES

For guidelines on cell sample preparation, please visit <https://www.aatbio.com/resources/guides/cell-sample-preparation.html>

SAMPLE EXPERIMENTAL PROTOCOL

Prepare and incubate cells with Apopxin™ Violet 450:

1. Treat cells with test compounds for a desired period of time (4 - 6 hours for Jurkat cells treated with staurosporine) to induce apoptosis.
2. Centrifuge the cells to get $1 - 5 \times 10^5$ cells/tube.
3. Resuspend cells in 200 µL of Assay Buffer (Component B).
4. Add 2 µL of Apopxin™ Violet 450 (Component A) into the cells.
5. **Optional:** Add 2 µL of 100X Propidium Iodide (Component C) into the cells for necrosis cells.
6. Incubate at room temperature for 30 to 60 minutes, protected from light.
7. Add 300 µL of Assay Buffer (Component B) to increase volume before analyzing the cells with a flow cytometer or fluorescence microscope.
8. Monitor the fluorescence intensity using a flow cytometer with 450/40 nm filter (Pacific Blue channel) or a fluorescence microscope with Violet filter.

Analyze by using a flow cytometer:

1. Quantify Apopxin™ Violet 450 binding using a flow cytometer with 450/40 nm filter (Pacific Blue channel). Measure the cell viability using 610/20 nm filter (PE-Texas Red channel) when propidium iodide is added into the cells.

Note Apopxin™ binding flow cytometric analysis on adherent cells is not routinely tested since specific membrane damage may occur during cell detachment or harvesting. However, methods for utilizing Annexin V for flow cytometry on adherent cell types have been previously reported by Casiola-Rosen et al. and van Engelen et al.

Analyze by using a fluorescence microscope:

1. Pipette the cell suspension after incubation, rinse 1 - 2 times with Assay Buffer, and then resuspend the cells with Assay Buffer.
2. Add the cells on a glass slide that is covered with a glass cover-slip.

Note For adherent cells, it is recommended to grow the cells directly on a cover-slip. After incubation with Apopxin™ Violet 450, rinse 1 - 2 times with Assay Buffer, and add Assay Buffer back to the cover-slip. Invert cover-slip on a glass slide and visualize the cells. The cells can also be fixed in 2% formaldehyde after the incubation with Apopxin™ Violet 450 and visualized under a microscope.

3. Analyze the apoptotic cells with Apopxin™ Violet 450 under a fluorescence microscope with Violet filter. Measure the cell viability using TRITC filter when propidium iodide is added into the cells. The blue staining on the plasma membrane indicates the Apopxin™ Violet 450 binding to PS on cell surface.

EXAMPLE DATA ANALYSIS AND FIGURES

In live non-apoptotic cells, Apopxin™ Violet 450 detects innate apoptosis in non-induced cells, which is typically 2-6% of all cells. In apoptotic cells Apopxin™ Violet

450 binds to phosphatidylserine, which is located on the outer leaflet of the cell membrane, resulted in increased staining intensity.




Figure 1. The detection of phosphatidylserine binding activity in Jurkat cells with Cell Meter™ Phosphatidylserine Apoptosis Assay Kit. Jurkat cells were treated without (Blue) or with 1 μM staurosporine (Red) in a 37 °C, 5% CO₂ incubator for 5 hours, and then dye loaded with Apopxin™ Violet 450 for 30 minutes. The fluorescence intensity of Apopxin™ Violet 450 was measured with a FACSCalibur (Becton Dickinson) flow cytometer using violet laser.

DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email info@aatbio.com if you have any questions.